

SIGAL MEDCO

Superior Quality Instruments



**ELECTROSURGICAL
INSTRUMENTS
CATALOGUE**

www.sigalmedco.com

For over 10 years SIGAL MEDCO has been devoted to providing its customers with the highest quality Electrosurgical, Surgical, Dental, Beauty and Podiatry Instruments.

SIGAL MEDCO is committed to provide you with the highest quality instruments. We are the original equipment manufacturer (OEM) for all above mention instruments.

Our production depth covers thousands of instruments ranging from the typical instruments use in clinic to sophisticated instruments required during surgeries. We deem, however, pertinent to submit further that the range and scope of our production activities is not limited to our production literature. We can produce instruments even if not depicted in our online catalogues. We receive samples, photos or references from the world-renowned companies catalogues.

QUALITY/ INNOVATION

At SIGAL MEDCO the quality is set in motion through motivated specifications planning by trained technical persons. The computer added design and imaging provide visual instruments a set configuration. This simplifies each aspect of instruments manufacturing.

QUALITY CERTIFICATIONS:

We are ISO 9001:2015, ISO 13485:2016, CE MARK, CGMP approved.

COMMITMENT:

Our commitment is to meet the changing needs of healthcare systems and to lead the industry by providing quality surgical instruments and dental instruments at competitive price.

We understand the diverse needs of our customers and know that quality along with delivery time are the most important and common requirements of our customers from both segments. This Realization allows us to serve our customers in the most efficient and cost effective manners.

Wish you all the Best.

SIGAL MEDCO Team

QUALITY CERTIFICATIONS

- ✓ ISO 9001:2008 Certified
- ✓ ISO 13485:2003 Certified
- ✓ CE Mark
- ✓ cGMP approved



ISO 9001:2015



ISO 13485:2016



CE Mark




cGMP



Surgical Instruments Association



Sialkot Chamber of Commerce



SIX SIGMA COMPANY

SIGAL MEDCO is a company believe in six-sigma mythology of quality control.

Our team consist of six sigma black belts who help control process variation.

FORGING TO PACKING

From forging to packing each step is controlled by our trained workers to make sure that each instruments manufactured is reflecting Hilmed's craftsmanship and commitments.

LATEST TECHNOLOGIES

Our production departments are well equipped with latest machines range from latest CNC machines, laser welding and NC machines to simple drilling machines to get the ultimate quality and process control.

European Bipolar Forceps	01 - 07
European Irrigation Bipolar Forceps	08
European Non-Stick Bipolar Forceps	09 - 12
USA 2 Pin Bipolar Forceps	13 - 16
UAA 2-Pin Non-Stick Bipolar Forceps	17 - 19
Silver Glide Non-Stick Bipolar Forceps	20 - 22
Single Use USA 2-Pin Bipolar Forceps	23 - 25
Single Use European Bipolar Forceps	26 - 27
Monopolar Cables	28
Bipolar Forceps Cables	29
Diathermy Instruments	30 - 31
Bipolar Artery Sealer	32
Bipolar Electrodes	33
Electrodes 2.4mm	34
Conization Electrodes	35
Electrodes 4.0mm	36
Patient Return Plate	37
Laparoscopic Electrodes	38
Electrosurgical Instruments for Gynecology	39 - 42

European Bipolar Forceps

McPherson (Str.)

L= 9.0 cm

Tip 0.5 mm



SM-8001

McPherson (crv.)

L= 11.5 cm

Tip 0.5 mm



SM-8002

Jeweler # 5

L= 11.5 cm

Tip 0.0 mm



SM-8003

Adson

L= 12.1 cm

Tip 0.1 mm



SM-8003

Jeweler # 4

L= 11.5 cm

Tip 0.5 mm



SM-8004

Jeweler # 7

L= 11.5 cm

Tip 0.5 mm



SM-8005

Coaptation

L= 10.8 cm

Tip 0.5 mm



SM-8006

Coaptation Nadler

L= 10.2 cm

Tip 0.5 mm



SM-8007

Iris (Str.)

L= 10.2 cm

Tip 0.5 mm



SM-8008

Iris (crv.)

L= 10.8 cm

Tip 1.0 mm



SM-8009

Straight

L= 16.0 cm

Tip 1.0 mm



SM-8010

Straight

L= 18.0 cm

Tip 1.0 mm



SM-8011

Straight

L= 20.0 cm

Tip 1.0 mm



SM-8012

Straight

L= 22.0 cm

Tip 1.0 mm



SM-8013

European Bipolar Forceps

Straight (Ang.)

L= 16.0 cm

Tip 1.0 mm



SM-8014

Straight (Ang.)

L= 18.0 cm

Tip 1.0 mm



SM-8015

Straight (Ang.)

L= 20.0 cm

Tip 1.0 mm



SM-8016

Straight (Ang.)

L= 22.0 cm

Tip 1.0 mm



SM-8017

Bayonet

L= 16.0 cm

Tip 1.0 mm



SM-8018

Bayonet

L= 18.0 cm

Tip 1.0 mm



SM-8019

Bayonet

L= 20.0 cm

Tip 1.0 mm



SM-8020

Bayonet (Ang.Up)

L= 16.0 cm

Tip 1.0 mm



SM-8021

Bayonet (Ang.Up)

L= 18.0 cm

Tip 1.0 mm



SM-8022

Bayonet (Ang.Up)

L= 20.0 cm

Tip 1.0 mm



SM-8023

Bayonet (Ang.Down)

L= 16.0 cm

Tip 1.0 mm



SM-8024

Bayonet (Ang.Up)

L= 18.0 cm

Tip 1.0 mm



SM-8025

Bayonet (Ang.Up)

L= 20.0 cm

Tip 1.0 mm



SM-8026

McPherson (Str.)

L= 9.0 cm

Tip 0.5 mm



SM-8027

European Bipolar Forceps

McPherson (crv.)

L= 11.5 cm

Tip 0.5 mm



SM-8028

Jeweler # 5

L= 11.5 cm

Tip 0.5 mm



SM-8029

Jeweler # 5

L= 11.5 cm

Tip 0.5 mm



SM-8030

Jeweler # 7

L= 11.5 cm

Tip 0.5 mm



SM-8031

Adson

L= 12.1 cm

Tip 1.0 mm



SM-8032

Straight

L= 16.0 cm

Tip 1.0 mm



SM-8033

Straight

L= 18.0 cm

Tip 1.0 mm



SM-8034

Straight

L= 20.0 cm

Tip 1.0 mm



SM-8035

Straight

L= 22.0 cm

Tip 1.0 mm



SM-8036

Straight (Ang.)

L= 16.0 cm

Tip 1.0 mm



SM-8037

Straight (Ang.)

L= 18.0 cm

Tip 1.0 mm



SM-8038

Straight (Ang.)

L= 20.0 cm

Tip 1.0 mm



SM-8039

Bayonet

L= 16.0 cm

Tip 1.0 mm



SM-8040

Bayonet

L= 18.0 cm

Tip 1.0 mm



SM-8041

European Bipolar Forceps

Bayonet

L= 22.0 cm

Tip 1.0 mm



SM-8042

Bayonet (Ang.Up)

L= 16.0 cm

Tip 1.0 mm



SM-8043

Bayonet (Ang.Up)

L= 18.0 cm

Tip 1.0 mm



SM-8044

Bayonet (Ang.Up)

L= 20.0 cm

Tip 1.0 mm



SM-8045

Bayonet (Ang.Down)

L= 16.0 cm

Tip 1.0 mm



SM-8046

Bayonet (Ang.Up)

L= 18.0 cm

Tip 1.0 mm



SM-8047

Bayonet (Ang.Up)

L= 20.0 cm

Tip 1.0 mm



SM-8048

Straight

L= 16.0 cm

Tip 1.0 mm



SM-8049 IR

Straight

L= 20.0 cm

Tip 1.0 mm



SM-8050 IR

Bayonet

L= 16.0 cm

Tip 1.0 mm



SM-8051 IR

Bayonet

L= 20.0 cm

Tip 1.0 mm



SM-8052 IR

European Non-Stick Bipolar Forcep

McPherson (Str.)

L= 9.0 cm

Tip 0.5 mm



SM-8053 NS

McPherson (crv.)

L= 11.5 cm

Tip 0.5 mm



SM-8054 NS

Jeweler # 5

L= 11.5 cm

Tip 0.5 mm



SM-8055 NS

Adson

L= 12.1 cm

Tip 1.0 mm



SM-8056 NS

Jeweler # 4

L= 11.5 cm

Tip 0.5 mm



SM-8057 NS

Jeweler # 7

L= 11.5 cm

Tip 0.5 mm



SM-8058 NS

Coaptation

L= 10.8 cm

Tip 0.5 mm



SM-8059 NS

Coaptation Nadler

L= 10.2 cm

Tip 0.5 mm



SM-8060 NS

Iris (Str.)

L= 10.8 cm

Tip 0.5 mm



SM-8061 NS

Iris (crv.)

L= 10.8 cm

Tip 1.0 mm



SM-8062 NS

Straight

L= 18.0 cm

Tip 1.0 mm



SM-8063 NS

Straight

L= 20.0 cm

Tip 1.0 mm



SM-8064 NS

Straight

L= 20.0 cm

Tip 1.0 mm




SM-8065 NS

Straight

L= 22.0 cm

Tip 1.0 mm



SM-8066 NS

European Non-Stick Bipolar Forcep

Straight (Ang.)

L= 16.0 cm

Tip 1.0 mm



SM-8067 NS

Straight (Ang.)

L= 18.5 cm

Tip 1.0 mm



SM-8068 NS

Straight (Ang.)

L= 20.0 cm

Tip 1.0 mm



SM-8069 NS

Straight (Ang.)

L= 22.0 cm

Tip 1.0 mm



SM-8070 NS

Bayonet

L= 16.0 cm

Tip 1.0 mm



SM-8071 NS

Bayonet

L= 18.0 cm

Tip 1.0 mm



SM-8072 NS

Bayonet

L= 20.0 cm

Tip 1.0 mm



SM-8073 NS

Bayonet (Ang.Up)

L= 16.0 cm

Tip 1.0 mm



SM-8074 NS

Bayonet (Ang.Up)

L= 18.0 cm

Tip 1.0 mm



SM-8075 NS

Bayonet (Ang.Up)

L= 20.0 cm

Tip 1.0 mm



SM-8076 NS

Bayonet (Ang.Down)

L= 16.0 cm

Tip 1.0 mm



SM-8077 NS

Bayonet (Ang.Up)

L= 18.0 cm

Tip 1.0 mm



SM-8078 NS

Bayonet (Ang.Up)








L= 20.0 cm

Tip 1.0 mm



SM-8079 NS

USA 2 Pin Bipolar Forceps

McPherson (Str.)		SM-8080
L= 9.0 cm		
Tip 0.5 mm		
McPherson (Crv.)		SM-8081
L= 9.0 cm		
Tip 0.5 mm		
Coaptation Nadler		SM-8082
L= 10.8 cm		
Tip 0.5 mm		
Coaptation		SM-8083
L= 10.8 cm		
Tip 0.5 mm		
Jeweler # 5		SM-8084
L= 11.0 cm		
Tip 0.5 mm		
Jeweler # 4		SM-8085
L= 11.5 cm		
Tip 0.5 mm		
Tenzel		SM-8087
L= 11.5 cm		
Tip 0.5 mm		

Jeweler # 7

L= 11.5 cm

Tip 0.5 mm



SM-8088

Iris (Str.)

L= 10.8 cm

Tip 0.5 mm



SM-8089

Iris (Crv.)

L= 10.8 cm

Tip 1.0 mm



SM-8090

Adson

L= 12.1 cm

Tip 1.0 mm



SM-8091

Semkin

L= 14.0 cm

Tip 1.0 mm



SM-8092

Cushing

L= 16.0 cm

Tip 1.0 mm



SM-8093

Cushing








L= 11.5 cm

Tip 1.0 mm



SM-8094

USA 2 Pin Bipolar Forceps

Cushing (Gursh)		SM-8095
L= 18.0 cm		
Tip 1.0 mm		
Gerald (Str)		SM-8096
L= 17.8 cm		
Tip 1.0 mm		
Gerald (Crv)		SM-8097
L= 17.8 cm		
Tip 1.0 mm		
Gerald Bayonet		SM-8098
L= 17.8 cm		
Tip 1.0 mm		
Jansen Bayonet		SM-8099
L= 15.0 cm		
Tip 1.0 mm		
Jansen Bayonet		SM-8100
L= 18.0 cm		
Tip 1.0 mm		
Malis		SM-8101
L= 17.8 cm		
Tip 1.0 mm		

Scollive Green-Wood

L= 19.1 cm

Tip 1.5 mm



SM-8102

Yasargil

L= 19.8 cm

Tip 1.0 mm



SM-8103

Hardy Bayonet

L= 12.5 cm

Tip 1.0 mm



SM-8104

Hardy Bayonet (Ang)








L= 21.5 cm

Tip 1.0 mm



SM-8105

USA 2-Pin Non-Stick Bipolar Forceps

McPherson (Str)		SM-8106 NS
L= 9.0 cm		
Tip 0.5 mm		
McPherson (Crv)		SM-8107 NS
L= 9.0 cm		
Tip 0.5 mm		
Jewelery #5		SM-8108 NS
L= 11.5 cm		
Tip 0.5 mm		
Jeweler # 4		SM-8109 NS
L= 11.5 cm		
Tip 0.5 mm		
Adson		SM-8110 NS
L= 12.1 cm		
Tip 1.0 mm		
Semkin		SM-8111 NS
L= 14.0 cm		
Tip 1.0 mm		
Cushing		SM-8112 NS
L= 16.0 cm		
Tip 1.0 mm		

Cushing

L= 18.0 cm

Tip 1.0 mm



SM-8113 NS

Gerald Bayonet

L= 17.8 cm

Tip 1.0 mm



SM-8114 NS

Jansen Bayonet

L= 18.0 cm

Tip 1.0 mm



SM-8115 NS

Scollive Green-Wood

L= 19.1 cm

Tip 1.5 mm



SM-8116 NS

Yasargil

L= 19.8 cm

Tip 1.0 mm



SM-8117 NS

Hardy Bayonet

L= 21.5 cm

Tip 1.0 mm



SM-8118 NS

USA 2 Pin Bipolar Forceps

Semkin

Irrigation



SM-8119 IR

Gerald Bayonet

Irrigation



SM-8120 IR

Jansen Bayonet

Irrigation



SM-8121 IR

Scollive Green-Wood

Irrigation



SM-8122 IR

Hardy Bayonet

Irrigation



SM-8123 IR

McPherson (Str.)

L= 9.0 cm

Tip 0.5 mm



SM-8124 FNS

McPherson (Crv.)

L= 9.0 cm

Tip 0.5 mm



SM-8125 FNS

Jeweler # 5

L= 11.5 cm

Tip 0.5 mm



SM-8126 FNS

Adson

L= 12.1 cm

Tip 1.0 mm



SM-8127 FNS

Cushing (Str.)

L= 16.0 cm

Tip 1.0 mm



SM-8128 FNS

Jansen Bayonet

L= 18.0 cm

Tip 1.0 mm



SM-8129 FNS

Jansen Bayonet

L= 20.0 cm

Tip 1.0 mm



SM-8130 FNS

Bayonet Full Non-Stick

L= 16.0 cm

Tip 1.0 mm



SM-8131 FNS

SilverGlide Non-Stick Bipolar Forcecep

Bayonet Full Non-Stick

L= 18.0 cm

Tip 1.0 mm



SM-8132 FNS

Bayonet Full Non-Stick

L= 20.0 cm

Tip 1.0 mm



SM-8133 FNS

Bayonet Full Non-Stick

L= 22.0 cm

Tip 1.0 mm



SM-8134 FNS

Bayonet (Irrigation)

L= 20.0 cm

Tip 1.0 mm



SM-8135 FNS

Bayonet (Irrigation)

L= 22.0 cm

Tip 1.0 mm



SM-8136 FNS

Bayonet Full Non-Stick

L= 22.0 cm

Tip 1.0 mm



SM-8137 FNS

Bayonet Full Non-Stick

L= 22.0 cm

Tip 0.5 mm



SM-8138 FNS

Bayonet Full Non-Stick

L= 24.0 cm

Tip 1.0 mm



SM-8139 FNS

Jeweler # 5

L= 11.5 cm

Tip 0.5 mm



SM-8140 FNS

Jeweler # 5

L= 11.5 cm

Tip 0.5 mm



SM-8142 FNS

Straight Full Non-Stick

L= 12.0 cm

Tip 1.0 mm



SM-8143 FNS

Straight (Ang.)

L= 16.0 cm

Tip 1.0 mm



SM-8144 FNS

Bayonet Full Non-Stick

L= 18.0 cm

Tip 1.0 mm



SM-8145 FNS

Bayonet Full Non-Stick

L= 20.0 cm

Tip 1.5 mm



SM-8146 FNS

Bayonet Full Non-Stick

L= 22.0 cm

Tip 1.5 mm



SM-8147 FNS

Bayonet Full Non-Stick

L= 22.0 cm

Tip 1.0 mm



SM-8148 FNS

Single Use USA 2-Pin Bipolar Forcep

Gerald Straight

Single Use



SM-8149 S

Jansen Bayonet

Irrigation



SM-8150 IRS

Jansen Bayonet

Single Use



SM-8151 S

Jeweler # 5

Single Use



SM-8152 S

Adson Forceps

Single Use



SM-8153 S

Cushing Forceps

Size: 16.0 cm

Tip: 1.0 mm



SM-8154 S

Jansen Bayonet ang-up

Size: 20.0 cm

Tip: 1.0 mm



SM-8154 S

Jeweler # 5

L= 11.5 cm

Tip 0.5 mm



SM-8155 FNS

Jeweler # 5

L= 11.5 cm

Tip 0.5 mm



SM-8156 FNS

Straight Full Non-Stick

L= 12.0 cm

Tip 1.0 mm



SM-8157 FNS

Straight (Ang.)

L= 16.0 cm

Tip 1.0 mm



SM-8158 FNS

Bayonet Full Non-Stick

L= 18.0 cm

Tip 1.0 mm



SM-8159 FNS

Bayonet Full Non-Stick

L= 20.0 cm

Tip 1.5 mm



SM-8160 FNS

Bayonet Full Non-Stick

L= 22.0 cm

Tip 1.5 mm



SM-8161 FNS

Bayonet Full Non-Stick

L= 22.0 cm

Tip 1.0 mm



SM-8162 FNS

Single Use USA 2-Pin Bipolar Forcep

Jansen Bayonet

Size: 18.0 cm

Tip: 1.0 mm



SM-8163 S

Jeweler # 4

Size: 11.5 cm

Tip: 0.5 mm



SM-8164 S

Adson Forceps

Size: 12.1 cm

Tip: 1.0 mm



SM-8165 S

Jansen Bayonet

Size: 18.0 cm

Tip: 1.0 mm



SM-8166 S

Jansen Bayonet

Size: 20.0 cm

Tip: 1.0 mm



SM-8167 S

Bayonet Forceps

Size: 18.0 cm

Tip: 1.0 mm



SM-8173 S

Bayonet Forceps

Size: 20.0 cm

Tip: 1.0 mm




SM-8174 S

Straging Forceps

Size: 20.0 cm

Tip: 1.0 mm



SM-8175 S

Straight Forceps (ang)

Size: 14.0 cm

Tip: 1.0 mm



SM-8176 S

Straight Forceps

Size: 14.0 cm

Tip: 1.0 mm



SM-8177 S

Single Use USA 2-Pin Bipolar Forceps

Jeweler # 5

Single Use



SM-8168 S

Straight Forceps

Single Use



SM-8169S

Bayonet Forceps

Single Use



SM-8170 S

Straight Forceps

Single Use



SM-8171 S

Straight Forceps (ang)

Single Use



SM-8172 S


Silicone Monopolar Cables 10" (3.0m)
SM-8185

Instrument Connector: Screw Type

Generator Connector: Standard Bovie 8.0mm Plug

Compatibility: Valleylab, Conmed, Bovie


Silicone Monopolar Cables 10" (3.0m)
SM-8186

Instrument Connector: 4.0mm Male Fitting

Generator Connector: Standard 4.0mm Banana Plug

Compatibility: Valleylab, Conmed & All ES Genrators


Silicone Monopolar Endoscoic Cable 10" (3.0m)
SM-8187

Instrument Connector: 4.0mm Female Plug

Generator Connector: Standard Bovie 8.0mm Plug

Compatibility: Valleylab, Conmed, Bovie


Silicone Monopolar Endoscopic Cable 10" (3.0m)
SM-8188

Instrument Connector: 4.0mm Female Plug

Generator Connector: Standard 4.0mm Banana Plug

Compatibility: Valleylab, Conmed & All ES Genrators


Silicone Diathermy Cable 10" (3.0m)
SM-8189

Instrument Connector: 4.8mm Male Fitting

Generator Connector: Standard Bovie 8.0mm Plug

Compatibility: Valleylab, Conmed, Bovie


Silicone Diathermy Cable 10" (3.0m)
SM-8190

Instrument Connector: 4.8mm Male Fitting

Generator Connector: Standard 4.0mm Banana Plug

Compatibility: Valleylab, Conmed & All ES Genrators

Bipolar Forceps Cables



Silicone Bipolar Cables 10" (3.0m)

SM-8178

Instrument Connector: USA 2-Pin Plug

Generator Connector: 4.0mm Banana Plug

Compatibility: Valleylab, Conmed, Codman, Olpympus and Mbst Bipolar Generator



Silicone Bipolar Cables 10" (3.0m)

SM-8179

Instrument Connector: European Flat Plug

Generator Connector: Martin Type

Compatibility: Aesculap, Martin, Wolf, Berchold



Silicone Bipolar Cables 10" (3.0m)

SM-8180

Instrument Connector: European Flat Plug

Generator Connector: 4.0mm Banana Plug

Compatibility: Valleylab, Conmed, Codman, Olpympus & Mbst Bipolar Generator



Silicone Bipolar Cables 10" (3.0m)

SM-8181

Instrument Connector: USA2-Pin Plug

Generator Connector: Martin Type

Compatibility: Aesculap, Martin, Wolf, Berchold



Silicone Bipolar Cables 10" (3.0m)

SM-8182

Instrument Connector: USA2-Pin Plug

Generator Connector: Erbe Type

Compatibility: Eibe, Siemens, Downs



Silicone Bipolar Cables 10" (3.0m)

SM-8183

Instrument Connector: USA2-Pin Plug

Generator Connector: Pin Distance 29.0mm

Compatibility: Valleylab, EMS, Dolley

Silicone Bipolar Cables 10" (3.0m)


SM-8184

Instrument Connector: European Flat Plug

Generator Connector: Pin Distance 29.0mm

Compatibility: Valleylab, EMS, Dolley






Turnur Warwick Scissor
(Straight)
SM-8191



Turnur Warwick Scissor
(Curved)
SM-8192



Surgical Scissor
(Straight)
SM-8193



Surgical Scissor
(Curved)
SM-8194



Abbey Sub Needle
17.5cm
SM-8195



Willson Hey Forceps
(Straight)
SM-8196

Richie Forceps
(Straight)
SM-8197



Willson Hey Forceps
(Curved)
SM-8198



Artery Forceps
(Straight)
SM-8199



Artery Forceps
(Curved)
SM-8200



Diathermy Instruments

Adson



SM-8201

Mcldone



SM-8202

Gillie



SM-8203

Dissecting Forceps



SM-8204

Gerald (1x2)



SM-8205

Adson Bayonet



SM-8206

Adson Bayonet

Angle Up



SM-8207

Remington Hob



SM-8208



Artery Forceps Bi-Clamp

SM-8209




Artery Forceps Non Stick

SM-8210




Artery Forceps

SM-8211



Artery Forceps

SM-8212



Silicone Bi-Clamp Connecting cord

Bipolar Electrodes



Straight

SM-8213

Angled

SM-8214

Angled

SM-8215

Angled 45°

SM-8216

Angled 90°

SM-8217



Bi-Turbinate Probe

Size: 17.5cm

SM-8218



Blade Electrode
Disposable: SM-8219



Ball Electrode (crv.) (3.0 mm)
Disposable: SM-8222



Needle Electrode
Disposable: SM-8225



Needle Electrode (crv.)
Disposable: SM-8228



Ball Electrode (2.0 mm)
Disposable: SM-8231



Ball Electrode (crv.) (2.0 mm)
Disposable: SM-8233



Blade Electrode (Shaft 13.0 cm)
Disposable: SM-8236



Blade Electrode (Shaft 13.0 cm)
Disposable: SM-8238



Blade Electrode (Shaft 13.0 cm)
Disposable: SM-8240



Ball Electrode (3.0 mm)
Disposable: SM-8220



Blade Electrode (crv.)
Disposable: SM-8223



Ball Electrode (4.0 mm)
Disposable: SM-8226



Ball Electrode (crv.) (4.0 mm)
Disposable: SM-8230



Ball Electrode (5.0 mm)
Disposable: SM-8232



Ball Electrode (crv.) (5.0 mm)
Disposable: SM-8234



Circle Loop Electrode 4.0 mm
Reusable: SM-8221



Circle Loop Electrode 4.0 mm
Reusable: SM-8224




Circle Loop Electrode 8.0 mm
Reusable: SM-8227



Diamond Electrode 8.0 mm
Reusable: SM-8229



Tip Cleaner SM-8235




Holster SM-8237



Blade Electrode (Shaft 13.0 cm)
Disposable: SM-8239


Conization Electrode
Disposable: SM-8241

Conization Electrodes




Conization Electrode

Disposable # SM-8242




Conization Electrode

Disposable # SM-8243



Conization Electrode

Disposable # SM-8244



Conization Electrode

Disposable # SM-8245



Conization Electrode

Disposable # SM-8246



Extender Electrode

SM-8247 Size: 7.0 cm



Extender Electrode

SM-8248 Size: 13.0 cm



Extender Electrode

SM-8249 Size: 18.0 cm



Blade Electrode
Reusable SM-8250



Blade Electrode (Crv)
Reusable SM-8252



Needle Electrode
Reusable SM-8254



Needle Electrode (Crv)
Reusable SM-8257



Ball Electrode (2mm)
Reusable SM-8259



Ball Electrode (Crv) (2mm)
Reusable SM-8262



Blade Electrode (Shaft 13.0cm)
Reusable SM-8265



Needle Electrode (Shaft 13.0cm)
Reusable SM-8267



Ball Electrode 3mm (Shaft 13.0 cm)
Reusable SM-8268



Ball Electrode 4mm (Shaft 13.0cm)
Reusable SM-8270



Ball Electrode 5mm (Shaft 13.0cm)
Reusable SM-8271



Ball Electrode
Reusable SM-8251



Ball Electrode (Crv) 3mm
Reusable SM-8253



Ball Electrode (4mm)
Reusable SM-8255



Ball Electrode (Crv) (4mm)
Reusable SM-8258



Ball Electrode (5mm)
Reusable SM-8260



Ball Electrode (Crv) (5mm)
Reusable SM-8263



Round Loop Electrode
(7.0cm) 4.0mm
Reusable SM-8256



Round Loop Electrode
(7.0cm) 6.0mm
Reusable SM-8259



Round Loop Electrode
(7.0cm) 8.0mm
Reusable SM-8261



Round Loop Electrode
(7.0cm) 10.0mm
Reusable SM-8264





Plate Electrode Small
Reusable SM-8266




Plate Electrode Large
Reusable SM-8269

Patient Return Plate

Patient Return Plate Reusable
5m Cable # SM-8271



Patient Return Plate Reusable
(5m Cable) SM-8272



Patient Return Cable
SM-8273- a(5m cable)
SM-8274- b(3m Cable)

Patient Return Cable
SM-8275- a(5m cable)
SM-8276- b(3m Cable)

S.S Patient Return Plate
(16x12cm) SM-8277

WITHOUT SUCTION


SM-8278-Hook Electrode (36cm)



SM-8280-Hook Electrode (36cm)



SM-8282-Needle-Hook Electrode (36cm)



SM-8284-Spatula-Hook Electrode (36cm)



SM-8286-Ball-Hook Electrode (36cm)

WITH SUCTION

 SM-8279L-Hook Electrode
(Irrigation) (36cm)

 SM-8281J-Hook Electrode
(Irrigation) (36cm)

 SM-8283-Needle-Hook Electrode
(Irrigation) (36cm)

 SM-8285-Spatula-Hook Electrode
(Irrigation) (36cm)

 SM-8287-Ball-Hook Electrode
(Irrigation) (36cm)

Laparoscopic Electrodes


Electrosurgical Instruments for Gynecology



Graves Vaginal Speculum
with smoke evacuation tube
Small: SM-8288
Medium: SM-8289
Large: SM-8290



VU-Max Speculum
with smoke evacuation tube
Small: SM-8291
Medium: SM-8292
Large: SM-8293



Pederson Vaginal Speculum
with smoke evacuation tube
Small: SM-8294
Medium: SM-8295
Large: SM-8296



VU-More (Graves) Speculum
with smoke evacuation tube
Small: SM-8300
Medium: SM-8301
Large: SM-8302

Cusco Speculum
with smoke evacuation tube
Small: SM-8297
Medium: SM-8298
Large: SM-8299



Collin Speculum
with smoke evacuation tube
Small: SM-8303
Medium: SM-8304
Large: SM-8305





**Weisman Speculum
with Suction (Left Open)**
Small: SM-8306
Medium: SM-8307
Large: SM-8308



**Weisman Speculum
with Suction (Right Open)**
Small: SM-8309
Medium: SM-8310
Large: SM-8311



Dressing Forceps
 SM-8312(21.0 cm)
 SM-8313(25.0 cm)




Tissue Forceps
 SM-8314(21.0 cm)
 SM-8315 (25.0 cm)



Iris Hook
 SM-8316(25.0 cm)



Emmet Hook
 SM-8317(25.0 cm)



Straight Hook
 SM-8318(25.0 cm)



Two-Prong Hook
 SM-8319(25.0 cm)



Electrosurgical Instruments for Gynecology



**Lateral Vaginal
Retractor - SM-8320**



**Kogan Endocervical
Speculum - SM-8321**



**Lateral Vaginal Retractor
(wide-Open) SM-8322**



**Rigby Retractor
SM-8323**



**Rigby Retractor
SM-8324**



3mm SM-8325



4mm SM-8326



3mm SM-8327



**Kogan Endocervical
Speculum**

SM-8329


**Lateral Vaginal
Retractor**

SM-8328



Graves Speculum
with smoking evacuation tube

**Small: SM-8330
Medium: SM-8331
Large: SM-8332**



Collin Speculum
with smoking evacuation tube

**Small: SM-8333
Medium: SM-8334
Large: SM-8335**



Pederson Graves Speculum
with smoking evacuation tube

**Small: SM-8336
Medium: SM-8337
Large: SM-8338**



Superior Quality Instruments

Galib Road Small Industrial Estate
Sialkot-51310 Pakistan

Tel/ Fax: 0092-52-3555075
What's app: 0092 332 8687141
What's app: 0092-300-2192008

Email: info@sigalmedco.com / sigalmedco@gmail.com
Web: <http://www.sigalmedco.com>
<http://www.sigaltools.com>
Alibaba Store: <http://sigalmedco.trustpass.alibaba.com>